

γ -eye

A highly sensitive,
benchtop, *in vivo*
imaging system for
all SPECT isotopes



**A novel imaging concept developed
specifically for accelerating
the early stages of preclinical
research.**

General

Our vision is to accelerate preclinical research, towards clinical translation for promising drugs, through our high-quality services and products.

We are a strong research partner who guides, consults, and supports all preclinical research studies of our collaborators.

- As a manufacturing company of novel breakthrough imaging systems (*eyes*), we offer simplicity, speed, and efficiency on a daily workflow, during the first steps of testing novel compounds.
- As a preclinical CRO, we offer a one-stop-shop at our state-of-the-art Laboratories that covers a full chain of preclinical studies, following a Good Laboratory Practice approach in the daily routine.

BIOEMTECH's ultimate goal is to aid scientists innovate, at every step of their research.



Introduction

BIOEMTECH's γ -eye™ is a novel imaging concept, developed specifically for accelerating the early stages of preclinical research. Suitable for a wide range of SPECT isotopes, γ -eye™ enables real time *in-vivo* imaging of mice from time zero post-injection while also it can provide accurately all temporal and spatial physiological distribution characteristics of the studied compounds. Designed based on the end-user's needs and integrated in an easy-to-use and intuitive environment, γ -eye™ is a unique imaging tool that combines high flexibility, efficiency, and accuracy for an extensive range of applications.

With a footprint of just 44 cm × 46 cm x 40 cm and a weight lower than 40 kg, γ -eye™ is a truly desktop device that can turn any space into an imaging lab. γ -eye™ comes with a laptop PC, which serves for data acquisition and processing. Standard licenses of the complete software suite Visual | eyes, are included.

Technology – Specifications



A. General Information

γ -eyeTM is a novel screening tool that offers the prosecution of fast, efficient and with high precision imaging studies. A large variety of SPECT isotopes as well as other imaginable signals resulting from alpha-emitting radionuclides can be studied in γ -eyeTM within seconds, without the cost of time-consuming post-processing routines. γ -eyeTM is a complete imaging solution, designed specifically to accelerate all stages of preclinical research.

Modality	SPECT
Anatomical mapping	Artificial X-ray
Active FOV	50 mm × 100 mm
Photodetectors	Photomultipliers Tubes
Scintillators	CsI:Na


B. Performance

γ -eyeTM, among other characteristics, infers good spatial resolution of 1.9 mm, energy resolution below 19%, while its high sensitivity enables real-time dynamic imaging with timeframes of 10 sec or lower (depending on the injection activity). Characterized by its high flexibility and based on the end user's needs, γ -eyeTM can be specifically optimized for certain applications and imaging studies.

Time frames	Down to 1 sec*
Sensitivity	341 cps/MBq
Spatial resolution	Up to 1.9 mm @ 0 mm
Energy resolution	Below 19%
Dynamic range	30 keV – 500 keV

C. Exchangeable collimators

γ -eyeTM, is designed to provide maximum flexibility to the user. Based on the type of the study, users can easily exchange collimators, choosing between a general purpose, a high-resolution, a high-sensitivity and a high-energy design. Specific designs can be provided upon request.




D. Anatomical Mapping

γ -eyeTM integrates an advance Artificial Intelligence Algorithm designed to synthesize morphological X-ray images by translating standard photographic images of mice. Artificially produced X-ray mouse images can be superimposed with functional radioisotope 2D images to enhance overall anatomical information.



Optical Photo



Real X-ray




Artificial X-ray

E. Animal Handling

To preserve animal's welfare and health, γ -eyeTM employs standard inputs for gas anaesthesia, fully compatible with third party systems. Anaesthesia then is provided into the mask of the animals, all throughout the imaging study. In addition, the system infers a heated imaging stage, thus maintaining the temperature of the animal at the desired level. Upon request, vital signs of the animal can be monitored including heart rate, respiratory rate, body temperature and oxygen saturation.



Anaesthesia	Standard inputs for gas anaesthesia; compatible with third party systems
Heating	Heated stage for optimum body temperature
Vital signs monitoring	Upon request



F. Footprint and connectivity

γ -eyeTM's footprint and standard digital interface connectivity can turn any space into an imaging lab. In addition, γ -eyeTM is characterized by simple power requirements and anaesthesia connections, allowing real time imaging and quantification inside a clean room, overcoming limitations, and facilitating scientists get great results under challenging conditions.

Outer dimensions	44 cm (L) × 46 cm (W) × 40 cm (H)
Weight	40 kg
AC input range	100-240 VAC
PC Connectivity	USB 2.0 Type A and GB Ethernet
Outer shielding	Sheet metal and acrylic




G. Visual Eyes Software


VISUAL | eyes is a complete software suite that serves for image acquisition, system control and analytical post-processing routines. Among other possibilities, users through VISUAL | eyes software, can generate imaging studies using custom and/or pre-defined protocols, obtain quantitative information in user's defined Region of Interests and export images to DICOM format.

Fast acquisition	Simplified procedure in a robust environment-
	Real time image visualization during the scan
Database	Raw data, DICOM storage, Compatibility with third party software
Anatomical mapping	Fusion with X-ray images artificially generated based on the mouse structural characteristics
Imaging protocols	Pre-defined and user's defined imaging protocols
Post processing	Integrated ROI manager for detailed post processing image analysis
License	Standard license for Mac and Windows

Live imaging console



Post Processing suite



Indicative Studies

- Tumor targeting for radiotherapy based on Pb212 radiopharmaceutical




Figure 1. 14 uCi of Pb212 – 20min imaging time

- Tumor targeting using I-125 radiopharmaceutical




Figure 2. 40 uCi of I-125 – 1h dynamic scan

- Inflammation imaging using Tc-99m compound





Figure 3. 20 uCi of Tc-99m paw injection – 20min imaging time



Scintica:

 info@scintica.com
www.scintica.com

London, ON

TEL: +1 (519) 914 5495
FAX: +1 (226) 884 5502

Webster, TX

TEL: +1 (832) 548 0895
FAX: +1 (226) 884 5502

London, England

TEL: +44 2045 252069
FAX: +1 (226) 884 5502


BIOEMTECH