



# $\phi$ -eye

A highly sensitive, benchtop,  
in vivo optical imaging  
system for preclinical studies



A highly sensitivity, benchtop, in vivo  
imaging system for preclinical optical  
studies

Scintica:

  
BIOEMTECH

# General

BIOEMTECH's vision is to accelerate preclinical research, towards clinical translation for promising drugs, through high quality services and products.

We are a strong research partner who guides, consults, and supports all preclinical research studies of our collaborators.

- As a manufacturing company of novel breakthrough imaging systems (*eyes*), we offer simplicity, speed, and efficiency on a daily workflow, during the first steps of testing novel compounds.
- As a preclinical CRO, we offer a one-stop-shop at our state-of-the-art Laboratories that covers a full chain of preclinical studies, following a Good Laboratory Practice approach in the daily routine.

BIOEMTECH's ultimate goal is to aid scientists innovate, at every step of their research.



## Introduction

BIOEMTECH's  $\phi$ -eye™ system is a high-sensitive, low-noise, *in vivo* optical imaging scanner. It enables non-invasive visualization and tracking of fluorescent and bioluminescent probes *in vivo*.  $\phi$ -eye™ is suitable for a wide range of applications, including but not limited to oncology. Cutting-edge technology combined with an intuitive, easy-to-use working environment, transforms  $\phi$ -eye™ into a unique imaging solution.

$\phi$ -eye™'s footprint of 60 cm × 60 cm x 60 cm and weight of lower than 60 kg, characterize it as a truly desktop device that can turn any place/space into an imaging lab. It comes with a laptop serving for data processing and acquisition while standard licenses of the complete software suite *Visual / eyes* are included.

## Technology- Specifications



## A. General Information

$\phi$ -eye<sup>TM</sup> enables non-invasive visualization and tracking of fluorescent and bioluminescent probes in the visible range within a living organism. A wide Field-of-View provides high throughput screening in a wide variety of preclinical studies. High expansion capability of fluorescence filter sets allows the conduction of studies with well-known and newly developed promising probes.

<b>Modality</b>	FLI / BLI
<b>Visible range</b>	400 – 800 nm
<b>Fluorescence filter sets</b>	4 (with high expansion capability)
<b>Number of mice</b>	Up to three (3)
<b>Active FOV</b>	17 cm × 17 cm
<b>White light</b>	White light for superior photographic images

## B. Performance


State-of-the-art technology offers superior sensitivity, resolution, wide Field-of-View and fast frame rates. Novel architecture and high-quality components provide superior fluorescence and bioluminescence images in a benchtop configuration. High quantum efficiency in the visible range enables the efficient screening of newly developed drugs *in vivo* using visible absorbing dyes.

<b>Detector</b>	sCMOS
<b>Resolution</b>	6.5-micron pixels, 2048 × 2048
<b>Quantum Efficiency</b>	82% @ 560 nm
<b>Readout noise (electrons)</b>	0.9 (median) / 1.5 (rms)
<b>Dark current (electrons/pixels/s)</b>	0.6
<b>Dynamic range (Typ.)</b>	33000:1
<b>Lens</b>	F/1.8
	>93% transmittance (VIS-NIR)
<b>Fluorescence bandpass filters</b>	>93% transmission
	>OD 6 blocking

### C. Animal handling

To preserve animals' welfare and health,  $\phi$ -eye<sup>TM</sup> employs standard inputs for gas anaesthesia - fully compatible with third party systems. Anaesthesia then is provided into the mask of the animals, all throughout the imaging study. In addition, the system infers a heated imaging stage, maintaining in such way the temperature of the animal in the desired level.



<b>Anaesthesia</b>	Standard inputs for gas anaesthesia; compatible with third party systems
<b>Heating</b>	Heated stage for optimum body temperature
<b>Light-tight imaging chamber</b>	30cm (L) × 30cm (W) × 30cm (H)



### D. Footprint and connectivity

$\phi$ -eye<sup>TM</sup>'s footprint and standard digital interface connectivity can turn any space into an imaging lab. In addition,  $\phi$ -eye<sup>TM</sup> is characterized by simple power requirements and anaesthesia connections, allowing real time imaging and quantification inside a clean room, overcoming limitations, and facilitating scientists get great results under challenging conditions.

<b>Outer dimensions</b>	60cm (L) × 60cm (W) × 60cm (H)
<b>Weight</b>	< 60 kg
<b>AC input range</b>	100-240 VAC
<b>PC Connectivity</b>	USB 3.0 and USB 2.0
<b>Outer shielding</b>	Sheet metal and acrylic




## E. Visual | eyes Software


The embedded software, *Visual | eyes* is a complete software suite, serving for image acquisition, system control and analytical post-processing routines. Among other possibilities, users through *Visual | eyes* software, can generate imaging studies using custom and/or pre-defined protocols, obtain quantitative information in user's defined Region of Interests and export images to DICOM format.

<b>Fast acquisition</b>	Simplified procedure in a robust environment
<b>Database</b>	Raw data, DICOM storage, Compatibility with third party software
<b>Dual layer imaging</b>	Superimposition with photographic image for superior anatomical mapping
<b>Imaging protocols</b>	Pre-defined and user's defined imaging protocols
<b>Post processing</b>	Integrated ROI manager for detailed post processing image analysis
<b>License</b>	Standalone licenses for Mac OS & Windows

# Live imaging console



# Post Processing suite




# Indicative Studies


## O Spondyloarthritis model

<b>Fluorescent dye</b>	Osteosense 680
<b>Mouse model</b>	TgA86 transgenic mouse
<b>Animal depilation</b>	Depilation of fur over hands, legs, and tail with cream, prior to imaging
<b>Administration route</b>	Retro-orbital
<b>Imaging</b>	23 hours post injection

### Imaging in prone position




M5\_M6\_M7 blank scan  
Exposure time 0,3sec




M5\_M6\_M7 Osteosense 680  
Exposure time 0,2sec




### Imaging in supine position



M5\_M6\_M7 blank scan  
Exposure time 0,3sec




M5\_M6\_M7 Osteosense 680  
Exposure time 0,2sec




### 0 Oncology study

<b>Fluorescent dye</b>	TdTomato
<b>Mouse model</b>	Genetically engineered mouse model (NF1-KO) that spontaneously develop cutaneous neurofibromas
<b>Animal depilation</b>	N/A
<b>Administration route</b>	N/A
<b>Imaging</b>	1-5 weeks



## ○ Cutaneous Neurofibromas model

Fluorescent dye	TdTomato
Mouse model	Genetically engineered mouse model (NF1-KO) that spontaneously develop cutaneous neurofibromas
Animal depilation	Depilation of fur over spine with cream, prior to imaging
Administration route	N/A
Imaging	1-5 weeks



## ○ Heart imaging after myocardial infarction


Fluorescent dye	Cy5
<b>Infarction surgery</b>	Thoracotomy and ligation at the level of the left-anterior descending (LAD) coronary artery (CA) approx. 4mm from the left auricle
<b>Animal depilation</b>	Depilation of fur over thorax and abdomen with cream, prior to imaging
<b>Administration route</b>	Intracardiac injection with 20uL (15 uM – Cy5)
<b>Imaging</b>	1 hour post injection

### Imaging in supine-prone-right lateral positions





**Scintica:**

 [info@scintica.com](mailto:info@scintica.com)  
[www.scintica.com](http://www.scintica.com)

**London, ON**

TEL: +1 (519) 914 5495  
FAX: +1 (226) 884 5502

**Webster, TX**

TEL: +1 (832) 548 0895  
FAX: +1 (226) 884 5502

**London, England**

TEL: +44 2045 252069  
FAX: +1 (226) 884 5502

  
BIOEMTECH